

MOBILE CLINIC with Full body DR check

Leading the nation's health
with mobile healthcare.



- **ASSURE:** Medical electrical appliances are inspected by Haier Biomedical
- **SPECIAL CUSTOMIZATION:** Customized Appearance Painting Service
- **FUNCTION:** support the inspection of multiple parts of the body
- **ACCEPTABLE CUSTOMIZATION:** equipment can be adjusted, purpose can be adjusted, design can be adjusted



扫描二维码关注海尔生物医疗
官方网站, 让生命更美好



扫描二维码关注海尔生物医疗
微信公众号, 让生命更美好

Qingdao Haier Biomedical Co., Ltd.
Hailemiao (Qingdao) Intelligent IoT Co., Ltd.
No. 280 Feng Yuan Road, High-tech Zone, Qingdao, 266109, P.R. China
URL: www.haierbiomedical.com
Global unified hotline: +86-532-88935593

Haier Biomedical
Make Life Better

Features

The vehicle design follows the basic principles of cabin optimization, electromechanical control, environmental control, and ergonomics, and is equipped with an imported energy-storage digital X-ray diagnosis and treatment system to meet the needs of all-position X-ray inspection.

9 meters square cabin structure, the protective wall adopts rigid steel skeleton structure, produced in strict accordance with the national lead room protection standard, 3.0 lead equivalent 6-sided protection, to ensure protection safety. , biochemical analysis and other inspections can be moved and popularized to the grassroots, leading the general health of the whole people with mobile medical care.



Parameters

Overall dimensions of the vehicle (mm)	11980*2500*3930
Vehicle weight (kg)	18995
Wheelbase	1800+5600mm
Chassis engine maximum net power (kW)	228
Maximum speed (km/h)	110



扫描二维码关注海尔生物医疗官方网站, 让生命更美好



扫描二维码关注海尔生物医疗微信公众号, 让生命更美好

Qingdao Haier Biomedical Co., Ltd.
 Hailemiao (Qingdao) Intelligent IoT Co., Ltd.
 No. 280 Feng Yuan Road, High-tech Zone, Qingdao, 266109, P.R. China
 URL: www.haierbiomedical.com
 Global unified hotline: +86-532-88935593

Haier Biomedical
 Make Life Better