

MOBILE LABORATORY

Mobile BSL-2S Protection Level Biosafety Laboratory



- **ASSURE:** Medical electrical appliances are inspected by Haier Biomedical
- **SPECIAL CUSTOMIZATION:** Customized Appearance Painting Service
- **AVOID INFECTION:** P2+ protection level
- **ACCEPTABLE CUSTOMIZATION:** equipment can be adjusted, purpose can be adjusted, design can be adjusted



扫描二维码关注海尔生物医疗
官方网站, 让生命更美好



扫描二维码关注海尔生物医疗
微信公众号, 让生命更美好

Qingdao Haier Biomedical Co., Ltd.
Hailemiao (Qingdao) Intelligent IoT Co., Ltd.
No. 280 Feng Yuan Road, High-tech Zone, Qingdao, 266109, P.R. China
URL: www.haierbiomedical.com
Global unified hotline: +86-532-88935593

Haier Biomedical
Make Life Better

Features

Haier Biomedical Mobile Biosafety Laboratory has the characteristics of flexibility, rapid response, safety and reliability, which can quickly reach around the epidemic area, and quickly carry out the collection, isolation and verification of pathogens. The car by the chassis and the combination with laboratory for negative pressure directional flow, can achieve BSL - 2S biological safety laboratory protection level, satisfy biosafety laboratory structure and facilities, equipment safety operation procedures and safety requirements, can satisfy the staff in handling contains pathogenic microorganism and toxins from the subjects to infect, surrounding environment from pollution.

transfer window	buffer chamber	biological safety cabinet
PCR machine	centrifuge	Chemical Analysis Instruments

Parameters

Overall dimensions of the vehicle (mm)	11980*2500*3960
Vehicle weight (kg)	17700
Wheelbase (mm)	1800*5600
Chassis engine maximum power (kW)	228
Maximum speed (km/h)	89



扫描二维码关注海尔生物医疗官方网站, 让生命更美好



扫描二维码关注海尔生物医疗微信公众号, 让生命更美好

Qingdao Haier Biomedical Co., Ltd.
Hailemiao (Qingdao) Intelligent IoT Co., Ltd.
No. 280 Feng Yuan Road, High-tech Zone, Qingdao, 266109, P.R. China
URL: www.haierbiomedical.com
Global unified hotline: +86-532-88935593

Haier Biomedical
Make Life Better