

Pharmacy Refrigerator

Innovative Design

- ◉ Purpose-built for medical and laboratory use
- ◉ Triple-paned glass to reduce condensation
- ◉ Eco-friendly hydrocarbon refrigerant and CFC-free insulation for reduced GWP
- ◉ Self-closing door hinge
- ◉ Multiple alarm functions for safer monitoring
- ◉ Ergonomic design with high performance, low energy consumption and low noise level



Low Energy Consumption
of 3kWh/day



Excellent Temperature
Uniformity of $\pm 1^{\circ}\text{C}$



Quiet Operation with
41dB(A)



HYC-1531GD

Qingdao Haier Biomedical Co., Ltd.

No.280 Fengyuan Road, High-tech Zone,
Qingdao, 266111, Shandong, P.R. China
E-mail: inquiry@haierbiomedical.com
Website: www.haiermedical.com

Oct, 2026



Haier Biomedical



Haier Biomedical
International



Haier Biomedical
International



Haier Biomedical
International



Haier Biomedical
International

Product Features



Two inverter compressors, dual cooling system ensure stable temperature control $\pm 1^{\circ}\text{C}$, low failure rate, and high energy efficiency 3kWh/24H.



Triple-paned tempered glass door with frame heater to reduce condensation even at 32°C , 80% RH ambient; self-closing door; magnetic gaskets to ensure a positive seal.



Multiple alarm functions for safer monitoring, audible/ visual alert with mute function; remote alarm contact is equipped as standard.



Anti-freezing thermostat is available to protect store items from freezing, which is DIN13277 compliance.



Energy efficient interior LED lighting, side-mounted to enlighten full space.



Two 25mm probe access ports on each side wall to accept user-provided monitoring equipment.



Optional USB module to easily record and export cabinet temperature and alarm data; optional RS485 port; optional printer to easily export cabinet temp data.



Adjustable plastic-coated wire shelves for improved circulation and flexible storage, 6 swivel casters with 2 leveling feet for easy movement and fixing.



Specifications

Model			HYC-1531GD
Technical Data	Cabinet Type		Upright, Glass Door
	Climate Class		N, 16~32°C
	Cooling Type		Forced Air Cooling
	Defrost Mode		Auto
	Sound Level (dB(A))		41
Performance	Temperature Range (°C)		2~8
	Set Temperature Range (°C)		4~15
	Temperature Uniformity (°C)		±1
Control	Controller		Microprocessor
	Display		LED
Electrical Data	Power Supply (V/Hz)		220~240/50/60; Dual Cooling System
Dimensions	Capacity (L/Cu.Ft)		1531/54.07
	Interior Dimensions (W*D*H)	mm	1750*594*1675
		in	68.90*23.39*65.94
	Exterior Dimensions (W*D*H)	mm	1840*730*1990
		in	72.44*28.74*78.35
	Packing Dimensions (W*D*H)	mm	1920*910*2110
		in	75.59*35.83*83.07
Alarms	High/Low Temperature		Y
	High Ambient Temperature		Y
	Door Ajar		Y
	Power Failure		Y
	Low Battery		Y
	Sensor Error		Y
	Communication Failure		Y
	Condenser Heat		Y
	Remote Alarm		Y
	Sound Light Alarm		Y
Accessories	Caster		6
	Foot		2
	Porthole		2
	Shelves/Drawers		18/-
	Printer		Optional
	USB Interface		Optional
	RS485		Optional
	Electromagnetic Lock		N/A
	Plastic Basket		Optional
Others	Certification		CE (Expected to obtain in 2026 Q2), DIN13277 (optional)

*Haier Biomedical reserves the right to change products and specifications without prior notice.